
# MAP SYSTEM

MICRO - APICAL PLACEMENT


Directions for use  
Mode d'emploi  
Gebrauchsanweisung  
указания по применению

Swiss quality  
dental products

problem

solved

---

**DIRECTIONS FOR USE****GB**

1. Indications for use.....	4
2. Contraindications .....	4
3. Warning.....	4
4. Precautions.....	4
5. Adverse reactions.....	4
6. Operating instructions.....	4
7. Disinfection, cleaning and sterilization.....	5

For dental professional use

## 1. Indications for use

The MAP system is used in hospital environments, clinics or dental offices by qualified dental personnel. The MAP System is designed for manually depositing any type of root canal repair material.

## 2. Contraindications

None known.

## 3. Warning

None known.

## 4. Precautions

4a. Sterilization of instruments before use is essential.

4b. Clean inside the head immediately before the repair material has hardened. If this is not done, it may be very difficult, if not impossible, to unblock and clean the device.

## 5. Adverse reactions

None known.

## 6. Operating instructions

**n** Refer to the illustrated user guide on the last page

### Preparing the instrument

- 1** Insert the piston into the corresponding head.  
Fix the selected head on the sleeve in the following way
- 2+3** Press the plunger so that the bayonet tip appears.
- 4** Fit the piston bayonet counter system and release the plunger.
- 5** Screw the head into the sleeve by twisting the body of the syringe.

### Filling

Ensure that at least one silicone stop is attached to the corresponding head (if the repair material needs to be applied successively, attach several silicone stops).  
Place the repair material in the well.

- 6** Thrust the end of the head into it several times so that it is filled (the plunger must be fully released so that the greatest possible volume is available).
- 7** Remove excess material by sliding the silicone stop outside the end of the head.

### Inserting the material

- 8** Press the plunger to expel the material.  
Fill the instrument again if necessary.
- 9** When the material has been fully inserted, press the plunger 2 or 3 times to expel as much material as possible.
- 10** Unscrew the head by twisting the body of the syringe.
- 11** Press the plunger of the syringe and take the piston\* out of the bayonet system.

\* PEEK pistons are intended for single use only

Clean the inside of the head using the cleaning curette. To clean the head thoroughly, insert the cleaning curette with a twisting, clockwise movement. Pull it out, again with a twisting, clockwise movement.

**12** Clean the parts with a cleaning and disinfecting product.

If necessary, clean the aluminum sterilization case (avoid using a dishwasher, ultrasonic bath and acid or basic chemicals for cleaning).

Dry the parts with compressed air.

Attach one or more silicone rings to the head (refills are available in dispensers, [REF](#) 20252 yellow, [REF](#) 20253 red, [REF](#) 20254 blue and [REF](#) 20255 green).

Place the parts in the sterilization case, without inserting the pistons in the corresponding heads.

Sterilize in the autoclave or by hot air.

## 7. Disinfection, cleaning and sterilization

Reprocessing procedure for dental instruments and radicular devices. Do not exceed 200°C (392°F)

### Foreword

For hygiene and sanitary safety purposes, all instruments not marked «sterile» must be cleaned, disinfected and sterilized before each usage to prevent any contamination.

### A. Area of application

Disinfection and Sterilization before first usage and reprocessing procedures concerning:

#### 1. Instruments

Root canal filling instruments (Pluggers, spreaders, compactors).

Supports, kits and instrument organizers.

#### B. Filling material

Please refer to the individual Direction for use.

### General recommendation

Use only a disinfecting solution which is approved for its efficacy (VAH/DGHM-listing, CE marking, FDA approval) and in accordance with the DFU of the disinfecting solution manufacturer. For all metal instruments, it is recommended to use anticorrosion disinfecting and cleaning agents.

For your own safety, please wear personal protective equipment (gloves, glasses, mask).

The user is responsible for the sterility of the product for the first cycle and each further usage as well as for the usage of damaged or dirty instruments where applicable after sterility.

Limitations and restrictions on reprocessing: The appearance of defects such as cracks, deformations (bent, twisted), corrosion, loss of color coding or marking, are indications that the devices are not able to fulfill the intended use with the required safety level.

Single use marked instruments are not approved for re-use.

The water quality has to comply with the local regulations especially for the last rinsing step or with a washer-disinfector.

Plastic supports, hand instruments and NiTi instruments are degraded by Hydrogen Peroxide (H<sub>2</sub>O<sub>2</sub>) solution.

NiTi Instruments are degraded if immersed more than 5 minutes in a solution of NaOCl at more than 5 %.

Instruments made of aluminum are degraded in presence of caustic soda solutions with mercury salt.

**Do not use** acid (pH < 6) or alkaline (pH > 8) solutions.

The washer-disinfector is not recommended for aluminum sterilization case.

## C. Step By Step Procedure

			A1. Instruments	
Operation	Operating mode	Warning	First use	Following uses
1	Disassembling	Disassemble the device, if required		✗
2	Pre-Disinfection	Soak all instruments immediately after use in a detergent and disinfecting solution combined with proteolytic enzyme if possible.	<ul style="list-style-type: none"> <li>• Silicone stops have to be removed.</li> <li>• Follow instructions and observe concentrations and immersion times given by the manufacturer (an excessive concentration may cause corrosion or other defects on instruments).</li> <li>• The disinfecting solution should be aldehyde free (to avoid blood impurities fixation) and without di- or triethanolamines as corrosion inhibitor.</li> <li>• Do not use disinfecting solutions containing Phenol or any products which are not compatible with the instruments (see general recommendations).</li> <li>• For visible impurities observed on instruments a pre-cleaning is recommended by brushing them manually with soft material.</li> </ul>	✗
3	Rinsing	Abundant rinsing (at least 1 min)	<ul style="list-style-type: none"> <li>• Use quality water in accordance with local regulations.</li> <li>• If a pre-disinfectant solution contains a corrosion inhibitor, it is recommended to rinse the instrument just before the cleaning.</li> </ul>	✗
4a	Automated cleaning with washer-disinfector	<ul style="list-style-type: none"> <li>• Place the devices in a kit, box or container to avoid any contact between instruments.</li> <li>• Put them in the washer-disinfector (Ao value &gt; 3000 or, at least 5 min at 90 °C).</li> </ul>	<ul style="list-style-type: none"> <li>• Discard any instruments with large obvious defects (broken, bent).</li> <li>• Avoid any contact between instruments. When placing in the washer-disinfector use kits, supports or container.</li> <li>• Follow instructions and observe concentrations given by the manufacturer (also see general recommendations).</li> <li>• Use only approved washer-disinfector according to EN ISO 15883, maintain and calibrate regularly.</li> </ul>	✗

			A1. Instruments		
	Operation	Operating mode	Warning	First use	Following uses
<b>OR</b>					
4b	Manual cleaning or assisted by an ultrasonic device	<ul style="list-style-type: none"> <li>Place the devices in a kit, box or container to avoid any contact between instruments.</li> <li>Immerse in the disinfecting solution with cleaning properties, assisted by an ultrasonic device if suitable.</li> </ul>	<ul style="list-style-type: none"> <li>No visible impurities should be observed on the instruments.</li> <li>Discard any instruments with large obvious defects (broken, bent, and twisted).</li> <li>Follow instructions and observe concentrations and time given by the manufacturer (see general recommendations under section 7).</li> <li>The disinfecting solution should be aldehyde free and without di- or triethanolamines as corrosion inhibitor.</li> </ul>	✗	✗
5	Rinsing	<ul style="list-style-type: none"> <li>Abundant rinsing (at least 1 min)</li> </ul>	<ul style="list-style-type: none"> <li>Use quality water in accordance with local regulations.</li> <li>If a disinfecting solution contains a corrosion inhibitor, it is recommended to rinse the instruments just before the autoclaving.</li> <li>Dry on a single use non-woven cloth, or with a drying machine or filtered compressed air.</li> </ul>	✗	✗
6	Inspection	<ul style="list-style-type: none"> <li>Inspect devices and sort out those with defects.</li> <li>Assemble the devices (stops)</li> </ul>	<ul style="list-style-type: none"> <li>Dirty instruments must be cleaned and disinfected again.</li> <li>Discard instruments which show any deformations (bent, twisted), damages (broken, corroded) or defects (loss of color coding or marking) affecting the resistance, the safety or the performance of the instrument.</li> </ul>	✗	✗
7	Packaging	<ul style="list-style-type: none"> <li>Place the devices in a kit, box or container to avoid any contact between instruments and pack the devices in "Sterilization pouches".</li> </ul>	<ul style="list-style-type: none"> <li>Check the validity period of the pouch given by the manufacturer to determine the shelf life.</li> <li>Use packaging which are resistant up to a temperature of 141°C (286°F) and in accordance with EN ISO 11607.</li> </ul>	✗	✗
8	Sterilization	<ul style="list-style-type: none"> <li>Steam sterilization at: 134 °C / 273°F during 18 min. at 2, 1 bar.</li> <li>Check the success of the sterilisable cycle (use physico-chemical indicator for each performed cycle)</li> </ul>	<ul style="list-style-type: none"> <li>The instruments, posts and the plastic supports must be sterilized according to the packaging labelling.</li> <li>Use only autoclaves that are matching the requirements of EN 13060, EN 285.</li> <li>The sterilization protocol has been validated by Produits Dentaires SA according to EN ISO 17665.</li> <li>Use only this recommended sterilization procedure.</li> <li>Control the efficiency (packaging integrity, no humidity, color change of sterilization indicators, physico-chemical integrators, digital records of cycles parameters).</li> <li>Traceability of procedure records.</li> </ul>	✗	✗
9	Storage	<ul style="list-style-type: none"> <li>Keep devices in sterilization packaging in a dry and clean environment</li> </ul>	<ul style="list-style-type: none"> <li>Sterility cannot be guaranteed if packaging is open, damaged or wet.</li> <li>Check the packaging and the medical devices before using them (packaging integrity, no humidity and validity period).</li> </ul>	✗	✗

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### Intro kit NiTi

REF 20294

For orthograde and retrograde obturations of anterior teeth

#### Contains

1 sterilisable aluminium case	185 x 95 x 30 mm
1 stainless syringe	
1 NiTi Memory Shape needle	N° 0 (yellow)
+ 3 spare PEEK plungers	N° 0
+ 1 cleaning curette	N° 0
1 curved needle (incl. NiTi plunger)	N° 0
1 NiTi Memory Shape needle	N° 1 (red)
+ 3 spare PEEK plungers	N° 1
+ 1 cleaning curette	N° 1
1 curved needle (incl. NiTi plunger)	N° 1
1 stainless steel well	ø 14 mm
1 paste cleaning device (key + stops)	N° 0
1 paste cleaning device (key + stops)	N° 1

### Universal kit NiTi

REF 20295

For orthograde and retrograde obturations of anterior and posterior teeth

#### Contains

1 sterilisable aluminium case	185 x 145 x 30 mm
1 stainless syringe	
1 NiTi Memory Shape needle	N° 0 (yellow)
1 Hooked needle	N° 0
+ 6 spare PEEK plungers	N° 0
+ 1 cleaning curette	N° 0
1 NiTi Memory Shape needle	N° 1 (red)
1 Hooked needle	N° 1
+ 6 spare PEEK plungers	N° 1
+ 1 cleaning curette	N° 1
1 curved needle (incl. NiTi plunger)	N° 0
1 curved needle (incl. NiTi plunger)	N° 1
1 stainless steel well	ø 14 mm
1 paste cleaning device (key + stops)	N° 0
1 paste cleaning device (key + stops)	N° 1

### Surgical kit

REF 20270

For retrograde obturations

#### Contains

1 sterilisable aluminium case	185 x 145 x 30 mm
1 stainless syringe	
1 triple-angle needle, right	N° 0/R (yellow)
1 triple-angle needle, left	N° 0/L
1 triple-angle needle, « retro » right	N° 0/R
1 triple-angle needle, « retro » left	N° 0/L
+ 6 spare PEEK plungers	N° 0
+ 1 cleaning curette	N° 0
1 triple-angle needle, right	N° 1/R (red)
1 triple-angle needle, left	N° 1/L
+ 6 pistons PEEK	N° 1
+ 1 cleaning curette	N° 1
1 stainless steel well	ø 14 mm
1 paste cleaning device (key + stops)	N° 0
1 paste cleaning device (key + stops)	N° 1

## PD MTA White

Endodontic filling material, specially developed to be placed with MAP System

### Advantages

- ▶ Optimized particle size
- ▶ Avoids bacterial migration
- ▶ Excellent marginal sealing capacity
- ▶ Stimulates the formation of a dentine layer (pulp capping)
- ▶ 2 sachets with 280 mg each

**Packaging**  
2 applications (2 × 280 mg)  
and 1 bottle of 3 ml pure H<sub>2</sub>O

REF  
20190

## Needles

### NITI «memory shape» needle

(2 PEEK plungers included)

N°	colour	ext. Ø mm	REF
0	yellow	0,90	20275
1	red	1,10	20277
2	blue	1,30	20278

### Straight needle (1 NITI plunger included)

N°	colour	ext. Ø mm	REF
0	yellow	0,90	20238
1	red	1,10	20239
2	blue	1,30	20292
3	green	1,50	20293

### Curved needle (1 NITI plunger included)

N°	colour	ext. Ø mm	REF
0	yellow	0,90	20236
1	red	1,10	20237
2	blue	1,30	20290
3	green	1,50	20291

### Triple-angle «retro» needle (1 PEEK plunger included)

N°	type	colour	ext. Ø mm	REF
0	right	yellow	0,90	20232
0	left	yellow	0,90	20233
1	right	red	1,10	20234
1	left	red	1,10	20235

### Triple-angle «retro» needle

(1 PEEK plunger included)

N°	type	colour	ext. Ø mm	REF
0	right	yellow	0,90	20271
0	left	yellow	0,90	20272
1	right	red	1,10	20273
1	left	red	1,10	20274

### Hooked needle (1 PEEK plunger included)

N°	colour	ext. Ø mm	REF
0	yellow	0,90	20263
1	red	1,10	20264
2	blue	1,30	20265
3	green	1,50	20266

## Accessories

### Plunger NITI

N°	colour	ext. Ø mm	REF
0	yellow	0,90	20247
1	red	1,10	20248
2	blue	1,30	20283
3	green	1,50	20284

### Plunger PEEK (6 pieces)

N°	colour	ext. Ø mm	REF
0	yellow	0,90	20245
1	red	1,10	20246
2	blue	1,30	20281
3	green	1,50	20282

### Cleaning curette NITI

N°	colour	ext. Ø mm	REF
0	yellow	0,90	20240
1	red	1,10	20241
2	blue	1,30	20243
3	green	1,50	20244

### Paste cleaning device

N°	colour	REF
0	yellow	20252
1	red	20253
2	blue	20254
3	green	20255