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
Polysiloxane, Siliciumdioxid, Pigmente, Additive, Platinkatalysator

Klassifikation

ISO 4823:2000, Type 2: Medium-bodied consistency

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
StatusBlue is an innovative A-silicone situation impression material. Short setting time in the mouth, easy removal from the mouth, unlimited shelf life and high detail reproduction make StatusBlue an alternative to alginates. StatusBlue is available in ready-to-use pre-filled MixStar- cartridges suitable for use in DMG MixStar and other commercially available mixing machines as well as in the 1:4 Automix safety-cartridge suitable for the Luxatemp applicator gun.

Indication

Impressions for provisionals, situation impressions, impressions of the opposing jaw, impressions for making vacuum-formed stents


Handling of the MixStar-Cartridge

Use in DMG MixStar-eMotion


1. Slide locking mechanism into open position. Do not remove locking mechanism from the cartridge.


2. Remove and discard the cap on the openings of base and catalyst outlets and attach a new mixing tip.

3. Secure the mixing tip by sliding the locking mechanism into locked position.

4. Insert the cartridge into the MixStar machine. A brief acoustical signal confirms the correct insertion of the cartridge.

- * Please be sure to follow the manufacturer's directions for the respective mixing device you are using.
- * Fresh cartridges are ready for use after the first 3 cm of impression material have been extruded and discarded.
- * After the application, leave the used tip on the cartridge to cap unused material.


Use in Other Mixing Units

If a MixStar cartridge is used with third party mixing devices insert the cartridge into the device first without the mixing tip. Otherwise appropriate mixing of the material cannot be guaranteed. In this case proceed as described in the section »Changing the mixing tip«, and follow the operating instructions for the device used.

Changing the mixing tip

Attach a new mixing tip before each mixing process. Remove the old mixing tip and any hardened material blocking the cartridge opening. First place the mixing tip on the hexagonal mixer shaft and secure it with a light turn, so that the shaft engages in the mixing tip. Then firmly attach the mixing.

Handling of the Automix-Cartridge

1. Insert cartridge

Push up the lever [1] at the back of the dispensing gun and pull back the slide [2] completely. Lift up the plastic lever [3], insert the cartridge and lock it by pushing down the lever [3]. Make sure that the notches on the cartridge line up with those on the dispensing gun.

2. Attach mixing tip

Remove the cartridge cap [4] or used tip by turning it counter-clockwise 90°, and discard it. Attach a new mixing tip [5]. Note the corresponding notches. Turn the mixing tip clockwise 90° to lock it.

3. Apply material

The material is mixed in the mixing tip by pulling the trigger [6] of the dispenser and can be applied directly.

When using a cartridge for the first time express a pea-size amount of material and discard it.

Leave the used tip on the cartridge to serve as a cap until the next use.

4. Remove cartridge

Push up the lever [1] at the back of the dispensing gun and pull back the slide [2] completely. Lift up the plastic lever [3] and remove the cartridge.

Recommended Use

Impression Tray

StatusBlue can be used in all common impression trays. To get better adhesion to metal or plastic trays, a tray adhesive may be applied.

Impression Taking

The prepared tray is loaded with the mixed material. To avoid air bubbles, the tip of the mixing tip should always stay immersed in the impression material. The tray is introduced into the mouth no later than 1:15 minutes after start of mixing. This position must be kept for at least 1:45 minutes. Remove the tray after setting.

Removal

Remove the impression parallel to the longitudinal axis of the prepared teeth. After removal, thoroughly clean with water and disinfect the impression.

Disinfection

Use only dipping solutions for the disinfection of the impression. Recommended is the use of a 2% glutaraldehyde solution. Please ensure that the used disinfectant is indicated for the use with addition curing silicone impression material.

After removal from the mouth, clean impression under running water for 15 seconds. Leave the impression in the disinfecting solution for about 10 to 15 minutes. Then rinse off the disinfection solution for 15 seconds under running water.

For more information, refer to the comments of the DGZMK. Comments of the DGZMK 3/93 V 1.0; www.dgzmk.de

ADA recommends following Council on Dental Materials, Instruments and Equipment. Disinfection of impressions. JADA 1991; 122(3): 110; www.ada.org

Casting and Galvanization

After disinfection, dry impression and store at normal room temperature (max. 78 °F / 25 °C). Impression should be poured 30 minutes at the earliest after removal from the mouth. We recommend using special purpose hard plaster (ISO 6873, Type 3; e.g. Moldano, Bayer) for the cast and an ultra hard dental plaster for the forms (ISO 6873, Type 4; e.g. Fujirock II, GC). StatusBlue impressions can be copper- or silverplated.

Storage of the Impression

Store impression at normal room temperature (78 °F / 25 °C). The impressions are stable and can be poured even after months without loss of accuracy.

Manufacture of Provisionals

Any provisional C&B material can be used. The setting of these materials is somewhat faster than with alginates.

Please Note

- Never allow the impression to come into contact with organic solvents or liquids containing same as this may result in swelling of the material and therefore may lead to inaccuracies.
- Possible residues of retraction aids can increase the setting time.
- Never allow the impression material to come into contact with latex gloves. Use of PVC or PE gloves is recommended.
- Do not mix StatusBlue with condensation curing silicone or polyether type impression materials.

Technical Data

StatusBlue-MixStar

Working time*	1:15 Min
Recommended time in mouth	1:45 Min
Dimensional change after 24 h	≤ 0.4%

Compression set	≤ 0.7%
Strain in compression	≈ 7.0%

StatusBlue-Automix

Working time*	1:15 Min
Recommended time in mouth	1:45 Min
Dimensional change after 24 h	≤ 0.4%
Compression set	≤ 0.7%
Strain in compression	≈ 7.0%

* Times refer to normal room temperature (72°F / 23°C) and a normal relative humidity of air (50%). Higher temperatures shorten, lower temperatures increase those times.

Composition

Addition curing polysiloxanes, silicon dioxide, food pigments, additives, platinum catalyst

Classification

ISO 4823:2000, Type 2: Medium-bodied consistency

Storage

Store and use in dry location at normal room temperature not exceeding 78 °F / 25 °C. We recommend not storing material in the refrigerator or in the freezer. Leave used mixing tip on the cartridge after use. It serves as a cap. Do not use after expiration date.

Packaging

StatusBlue-MixStar

- 1 MixStar-cartridge @ 380 ml paste
- 10 MixStar-Tips REF 909693
USA, please use REF 999693
- 5 MixStar-cartridges @ 380 ml paste
- 50 MixStar-Tips REF 909694

StatusBlue-Automix

- 2 Automix-cartridges @ 50 ml paste
- 4 Automix-Tips REF 909790

- 8 Automix-cartridges @ 50 ml paste
- 16 Automix-Tips REF 909791
USA, please use REF 999826

Caution: Federal Law restricts this device to sale by, or on the order of a dentist, or other practitioner licensed by law of the state in which he or she practices to use or order the use of this device.

Keep away from children!
For dental use only!

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